

Product Datasheet

TRYPTOPHAN HYDROXYLASE antibody (orb431996)

Description	Sheep polyclonal antibody to TRYPTOPHAN HYDROXYLASE
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, WB
Immunogen	Recombinant rabbit tryptophan hydroxylase purified from E. coli.
Target	TRYPTOPHAN HYDROXYLASE
Preservatives	0.09% Sodium Azide 0.01% Bovine Serum Albumin 50% Glycerol. 10mM HEPES pH 7.5
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.