

## Product Datasheet

### CD290 antibody (orb432001)

|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CD290                                             |
| <b>Species/Host</b>        | Rabbit                                                                          |
| <b>Reactivity</b>          | Human                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | WB                                                                              |
| <b>Immunogen</b>           | Synthetic peptide corresponding to a 14 amino acid sequence within human CD290. |
| <b>Target</b>              | CD290                                                                           |
| <b>Preservatives</b>       | 0.02% Sodium Azide. Phosphate buffered saline                                   |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                           |
| <b>Concentration</b>       | 1 mg/ml                                                                         |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                        |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | Polyclonal IgG                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                      |
| <b>Purity</b>              | Purified                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |