

Product Datasheet

IFN ALPHA antibody (orb432010)

Description	Rabbit polyclonal antibody to IFN ALPHA
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, FA, IP, WB
Target	IFN ALPHA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.