

Product Datasheet

GREEN FLUORESCENT PROTEIN antibody (orb432030)

Description	Goat polyclonal antibody to GREEN FLUORESCENT PROTEIN
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Fusion protein corresponding to full length GFP derived from Aequorea victoria
Target	GREEN FLUORESCENT PROTEIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.