

Product Datasheet

GAPDH antibody (orb432042)

Description	Rabbit polyclonal antibody to GAPDH
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A 16 amino acid peptide located near human GAPDH carboxy-terminus.
Target	GAPDH
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.