

Product Datasheet

ST2 antibody (orb432045)

Description	Rabbit polyclonal antibody to ST2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A 16 amino acid synthetic peptide located near mouse ST2 amino-terminus. This peptide is common to the three known isoforms of ST2.
Target	ST2
Preservatives	0.02% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.