

Product Datasheet

CD4/CD8 antibody (FITC/RPE) (orb432067)

Description	Mouse monoclonal antibody against CD4/CD8 conjugated to FITC/RPE
Species/Host	Rat
Reactivity	Canine
Conjugation	FITC/RPE
Tested Applications	FC
Target	CD4/CD8
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Dual Colour combination consisting of FITC conjugated and RPE conjugated monoclonal antibodies mixed in optimal ratio - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Cocktail
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.