

Product Datasheet

CD44 antibody (orb432078)

Description	Mouse monoclonal antibody to CD44
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Spleen cells from B6PL-Ly-2a / Ly-3a B6-Ly-1a mice.
Target	CD44
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.