

Product Datasheet

CD4 antibody (orb432087)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC
Target	CD4
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(DOMAINS 3 AND 4)
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.