

Product Datasheet

OX-62 antibody (orb432095)

Description	Mouse monoclonal antibody to OX-62
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP, WB
Immunogen	Density gradient enriched PVG rat thoracic duct dendritic cells.
Target	OX-62
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.