

Product Datasheet

CD5 antibody (APC) (orb432114)

Description	Rat monoclonal antibody against CD5 conjugated to APC
Species/Host	Rat
Reactivity	Canine
Conjugation	APC
Tested Applications	FC
Immunogen	Concanavilin A activated canine peripheral blood cells
Target	CD5
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Allophycocyanin (APC) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.