

Product Datasheet

MHC CLASS II MONOMORPHIC antibody (FITC) (orb432133)

Description	Rat monoclonal antibody against MHC CLASS II MONOMORPHIC conjugated to FITC
Species/Host	Rat
Reactivity	Canine
Conjugation	FITC
Tested Applications	FC
Immunogen	Con A activated canine peripheral blood cells.
Target	MHC CLASS II MONOMORPHIC
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.