

Product Datasheet

MHC CLASS II MONOMORPHIC antibody (orb432134)

Description	Rat monoclonal antibody to MHC CLASS II MONOMORPHIC
Species/Host	Rat
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	Con A activated canine peripheral blood cells.
Target	MHC CLASS II MONOMORPHIC
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.