

Product Datasheet

CD59 antibody (orb432148)

Description	Mouse monoclonal antibody to CD59
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, IHC-P, IP, WB
Immunogen	Thymocytes and T lymphocytes.
Target	CD59
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.