

Product Datasheet

CD32 antibody (orb432150)

Description	Mouse monoclonal antibody to CD32
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr, IP
Immunogen	K562 cell line.
Target	CD32
Preservatives	0.09% sodium azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.