

Product Datasheet

TCR GAMMA/DELTA antibody (orb432234)

Description	Mouse monoclonal antibody to TCR GAMMA/DELTA
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	CD3 +ve TCR alpha/beta -ve rat T cell hybridoma line.
Target	TCR GAMMA/DELTA
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.