

Product Datasheet

CD71 antibody (orb432239)

Description	Mouse monoclonal antibody to CD71
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF
Immunogen	KGI cell line.
Target	CD71
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.