

Product Datasheet

PI-3 KINASE p85 SUBUNIT ALPHA antibody (orb432246)

Description	Mouse monoclonal antibody to PI-3 KINASE p85 SUBUNIT ALPHA
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	IH-Fr, IP, WB
Immunogen	Baculovirus expressed (Recombinant) Bovine p85 alpha subunit of PI 3 kinase.
Target	PI-3 KINASE p85 SUBUNIT ALPHA
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.