

## Product Datasheet

### PI-3 KINASE p85 SUBUNIT BETA antibody (orb432249)

|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to PI-3 KINASE p85 SUBUNIT BETA                    |
| <b>Species/Host</b>        | Mouse                                                                        |
| <b>Reactivity</b>          | Bovine                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                 |
| <b>Tested Applications</b> | IH-Fr, IP, WB                                                                |
| <b>Immunogen</b>           | Balculovirus expressed (Recombinant) p85 beta subunit of bovine PI-3 kinase. |
| <b>Target</b>              | PI-3 KINASE p85 SUBUNIT BETA                                                 |
| <b>Preservatives</b>       | 0.09% sodium azide (NaN <sub>3</sub> ).<br>Phosphate buffered saline         |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                        |
| <b>Concentration</b>       | 1.0 mg/ml                                                                    |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                     |
| <b>Note</b>                | For research use only                                                        |
| <b>Isotype</b>             | IgG1                                                                         |
| <b>Clonality</b>           | Monoclonal                                                                   |
| <b>Purity</b>              | Purified                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                              |