

Product Datasheet

CD7 antibody (orb432272)

Description	Mouse monoclonal antibody to CD7
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	Normal human peripheral blood lymphocytes.
Target	CD7
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.