

Product Datasheet

TROPONIN I antibody (orb432277)

Description	Mouse monoclonal antibody to TROPONIN I
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human Cardiac troponin I.
Target	TROPONIN I
Preservatives	0.1% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(CARDIAC)
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.