

Product Datasheet

SWC8 antibody (orb432280)

Description	Mouse monoclonal antibody to SWC8
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr
Immunogen	Porcine Lamina Propria Leucocytes.
Target	SWC8
Preservatives	0.09% sodium azide (NaN ₃)
Form/Appearance	Tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.