

Product Datasheet

CD11R1 antibody (orb432281)

Description	Mouse monoclonal antibody to CD11R1
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	Porcine Lamina Propria Leucocytes.
Target	CD11R1
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.