

## Product Datasheet

### CD45R antibody (orb432313)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Rat monoclonal antibody to CD45R                                  |
| <b>Species/Host</b>        | Rat                                                               |
| <b>Reactivity</b>          | Mouse                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | FACS, ICC, IF, IHC-Fr, IHC-P, IP                                  |
| <b>Immunogen</b>           | Murine leukemia-induced pre-B tumor cells (RAW112)                |
| <b>Target</b>              | CD45R                                                             |
| <b>Preservatives</b>       | 0.09% sodium azide (NaN <sub>3</sub> ). Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                             |
| <b>Concentration</b>       | 1.0mg/ml                                                          |
| <b>Storage</b>             | +4°C, -20°C if preferred                                          |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG2a                                                             |
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>Purity</b>              | Purified                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |