

Product Datasheet

CD2 antibody (orb432365)

Description	Rat monoclonal antibody to CD2
Species/Host	Rat
Reactivity	Equine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Equine peripheral blood lymphocytes.
Target	CD2
Preservatives	0.09% sodium azide (NaN ₃)
Form/Appearance	Tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.