

Product Datasheet

THYROID STIMULATING HORMONE RECEPTOR antibody (orb432366)

Description	Mouse monoclonal antibody to THYROID STIMULATING HORMONE RECEPTOR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP, WB
Immunogen	Recombinant Human TSH receptor.
Target	THYROID STIMULATING HORMONE RECEPTOR
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.