

Product Datasheet

CD204 antibody (orb432397)

Description	Rat monoclonal antibody to CD204
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, IP, WB
Immunogen	Raw 264 cell line
Target	CD204
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.