

Product Datasheet

FITC antibody (orb432411)

Description	Rat monoclonal antibody to FITC
Species/Host	Rat
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	FITC conjugated mouse immunoglobulin.
Target	FITC
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgM - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.