

Product Datasheet

CD31 antibody (orb432414)

Description	Mouse monoclonal antibody to CD31
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, WB
Immunogen	Activated, Lewis rat derived microglial cells.
Target	CD31
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.