

Product Datasheet

CD8 ALPHA/BETA antibody (orb432433)

Description	Mouse monoclonal antibody to CD8 ALPHA/BETA
Species/Host	Mouse
Reactivity	Feline
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Feline thymocytes.
Target	CD8 ALPHA/BETA
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.