

Product Datasheet

V5-TAG antibody (orb432443)

Description	Mouse monoclonal antibody to V5-TAG
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, IP, WB
Immunogen	Paramyxovirus Simian-Virus 5 (SV5)
Target	V5-TAG
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.