

Product Datasheet

CD117 antibody (orb432446)

Description	Rat monoclonal antibody to CD117
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	IL-3 dependent mast cells.
Target	CD117
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.