

Product Datasheet

CD11c antibody (orb432453)

Description	Hamster monoclonal antibody to CD11c
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IP
Immunogen	Mouse spleen dendritic cells.
Target	CD11c
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.