

Product Datasheet

CD11c antibody (Biotin) (orb432454)

Description	Hamster monoclonal antibody against CD11c conjugated to Biotin
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Mouse spleen dendritic cells.
Target	CD11c
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.