

Product Datasheet

CD11c antibody (RPE) (orb432458)

Description	Hamster monoclonal antibody against CD11c conjugated to RPE
Species/Host	Hamster
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Mouse spleen dendritic cells.
Target	CD11c
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin5% sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.