

Product Datasheet

CD31 antibody (orb432461)

Description	Hamster monoclonal antibody to CD31
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FACS, ICC, IF
Immunogen	D10.G4.1 cells (Kaye et al. 1984).
Target	CD31
Preservatives	None present.Sterile filtered.. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.