

Product Datasheet

HISTIDINE TAG antibody (Biotin) (orb432479)

Description	Mouse monoclonal antibody against HISTIDINE TAG conjugated to Biotin
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	PAX6 transcription factor linked to histidine tag.
Target	HISTIDINE TAG
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.