

Product Datasheet

FIBROBLASTS/EPITHELIAL CELLS antibody (orb432485)

Description	Mouse monoclonal antibody to FIBROBLASTS/EPITHELIAL CELLS
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	Human foreskin fibroblasts.
Target	FIBROBLASTS/EPITHELIAL CELLS
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.