

Product Datasheet

TCR ALPHA/BETA antibody (FITC) (orb432493)

Description	Hamster monoclonal antibody against TCR ALPHA/BETA conjugated to FITC
Species/Host	Hamster
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Affinity purified murine T cell receptor.
Target	TCR ALPHA/BETA
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	BETA CHAIN
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.