

Product Datasheet

CD161 antibody (FITC) (orb432514)

Description	Mouse monoclonal antibody against CD161 conjugated to FITC
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Purified splenic NK cells from the LEW rat strain.
Target	CD161
Preservatives	0.09% sodium azide (NaN ₃)1% bovine serum albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.