

Product Datasheet

CD161 antibody (orb432519)

Description	Mouse monoclonal antibody to CD161
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr, IP, WB
Immunogen	Purified splenic NK cells from the LEW rat strain.
Target	CD161
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.