

Product Datasheet

MAST CELL TRYPTASE antibody (orb432524)

Description	Mouse monoclonal antibody to MAST CELL TRYPTASE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Human mast cell tryptase purified from human lung tissue.
Target	MAST CELL TRYPTASE
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.