

Product Datasheet

AGGRECAN antibody (orb432543)

Description	Mouse monoclonal antibody to AGGRECAN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, IHC-P, IP, WB
Immunogen	Purified human articular cartilage aggrecan
Target	AGGRECAN
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.