

Product Datasheet

CARTILAGE OLIGOMERIC MATRIX PROTEIN antibody (orb432544)

Description	Rat monoclonal antibody to CARTILAGE OLIGOMERIC MATRIX PROTEIN
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, IHC-P, IP, WB
Immunogen	Human cartilage derived COMP
Target	CARTILAGE OLIGOMERIC MATRIX PROTEIN
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.