

Product Datasheet

TNF ALPHA antibody (orb432573)

Description	Rat monoclonal antibody to TNF ALPHA
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr
Immunogen	Recombinant murine TNF alpha expressed in E. coli.
Target	TNF ALPHA
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.