

Product Datasheet

IL-4 antibody (orb432588)

Description	Rat monoclonal antibody to IL-4
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, WB
Immunogen	Purified CHO cell expressed recombinant human IL-4.
Target	IL-4
Preservatives	0.09% sodium azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.