

Product Datasheet

KAPPA LIGHT CHAIN antibody (HRP) (orb432589)

Description	Rat monoclonal antibody against KAPPA LIGHT CHAIN conjugated to HRP
Species/Host	Rat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Purified mouse IgG.
Target	KAPPA LIGHT CHAIN
Preservatives	0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.