

Product Datasheet

CD95 antibody (orb432594)

Description	Mouse monoclonal antibody to CD95
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Fusion protein comprising extracellular domain of human Fas linked to human Fc.
Target	CD95
Preservatives	0.1% sodium azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.