

Product Datasheet

CD4 antibody (orb432599)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	MLR generated rat T helper lymphocytes.
Target	CD4
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(DOMAIN 2)
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.