

Product Datasheet

Bcl-2 antibody (orb432609)

Description	Mouse monoclonal antibody to Bcl-2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, IP, WB
Immunogen	Synthetic peptide, amino acids 41-54 of the Bcl-2 protein
Target	Bcl-2
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.